

whiz[®] Relief bag


PRODUCT LAUNCH:

The Relief bag is a urine collection bag for use with the Whiz[®] freedom. It is ideal for emergencies, traffic jams, medical needs, or simply when it's too cold to go outside!

Easy to Use: (See reverse of bag for detailed instructions):

1. Open the clip, squeeze the narrow end of the Whiz[®] together and push it 1 to 2 cm into the Relief Bag tube. It will be a tight fit.
2. Position the Relief Bag below the Whiz[®] and wee normally through the Whiz[®].
3. When finished firmly close the clip and remove the Whiz[®]
4. Dispose of urine at your first opportunity after use.

The Relief Bag is made from 250 micron thickness PVC and a Silicon tube. BS EN ISO 8669-2:1997 Compliant.
Capacity: 1.6L.



Whiz[®] Relief Bag



RRP: £ 3.95

Combo Pack: freedom & Relief Bag



£ 12.50

Ultimate Pack: Whiz[®] freedom + Whiz Bag + Relief Bag, 10cm Tubing + 15ml Hand Sanitizer



£ 19.95

The Relief Bag can be washed out and re-used. However, it is primarily designed for emergency and one-off situations and should not be re-used on a regular basis over an extended period of time.

ARTWORK: Images of all Whiz[®] freedom product, accessories and sets can be downloaded from: <ftp://213.171.193.5/> username: whizart password: whiz2009

Further Information & Ordering:

Please contact Guy or Amaury
Amaury: a.marinho@jbol.co.uk

Whiz[®] is a registered trademark of Jbol Ltd.

Jbol Ltd t/a Whiz UK
1 Folly Bridge
OX14LB
Tel: +44 (0) 1865 240572
Fax: +44 (0) 1865 246565

Guy Burnell
Business Development Manager
guy@jbol.co.uk
Mob: + 44 (0) 7770 333395
www.whizfreedom.com